


FTIR/FT-NIR bench top spectrometer Interspec 200-X



Specifications:

Wavelength range, IR	7000 to 400 cm^{-1}
Wavelength range, NIR	15000 to 3850 cm^{-1}
Resolution, standard	2 cm^{-1} (8 cm^{-1} - NIR)
Resolution, option	1.0 or 0.5 cm^{-1}
Sample compartment	W20xD26xH16 cm
Beam at sample	10 mm dia.
Operating system	Windows 98/XP
Power	90-230 VAC, 12 VDC, 40 W
Dimensions	W59xD39xH19 cm
Weight	24 kg

Interspectrum OU - Toravere, 61602 Tartumaa, Estonia
Tel.: +372 7383008, Fax: +372 7383065
Email: info@interspectrum.ee,
www.interspectrum.ee