

Portable FTIR/FT-NIR spectrometer Interspec 300-X



Specifications:

Wavelength range, IR	7000 to 400 cm^{-1} (without ATR accessory)
Wavelength range, NIR	15000 to 3850 cm^{-1}
Resolution, standard	2 cm^{-1} (8 cm^{-1} - NIR)
Resolution, option	1.0 or 0.5 cm^{-1}
Operating system	Windows 98/XP
Power	90-230 VAC, 12 VDC, 40 W
Dimensions	W49xD39xH20 cm
Weight	18 kg
Temp. environment	15 – 28 $^{\circ}\text{C}$
Humidity environment	Best below 65%

Interspectrum OU - Toravere, 61602 Tartumaa, Estonia
Tel.: +372 7383008, Fax: +372 7383065
Email: info@interspectrum.ee,
www.interspectrum.ee