

FTIR/FT-NIR module Interspec 400-X



Specifications:

Wavelength range, IR	7000 to 400 cm ⁻¹
Wavelength range, NIR	15000 to 3000 cm ⁻¹
Resolution, standard	2 cm ⁻¹ (8 cm ⁻¹ - NIR)
Resolution, option	1.0 or 0.5 cm ⁻¹
Output beam diameter	38 mm
Interface	Arcnet
Power	2.35 VDC (5.5A), +/-12 V, +5 V, 35 W
Dimensions	W33xD25xH16 cm
Weight	11 kg

Interspectrum OU - Toravere, 61602 Tartumaa, Estonia
Tel.: +372 7383008, Fax: +372 7383065
Email: info@interspectrum.ee
www.interspectrum.ee