

FTIR/FT-NIR module Interspec 401-X



Specifications:

Wavelength range, IR	7000 to 400 cm^{-1}
Wavelength range, NIR	15000 to 3000 cm^{-1}
Resolution, standard	2 cm^{-1} (8 cm^{-1} - NIR)
Resolution, option	1.0 or 0.5 cm^{-1}
Output beam diameter	38 mm
Data acquisition system	16 bit
Power	12 V
Interface	Arcnet
Dimensions	W38xD25xH16 cm
Weight	12 kg

Interspectrum OU - Toravere, 61602 Tartumaa, Estonia
Tel.: +372 7383008, Fax: +372 7383065
Email: info@interspectrum.ee,
www.interspectrum.ee